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DIRECTIONS**Light Curing material****1. General terms:****1.1. Name of the product: Light Curing Material**

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1.4. Storage conditions:

- cool, protected from light and frost
- min. -10 °C to (medium term lasting) max.+ 50°C

2. Determined use:

Synthetic material easy to adapt at room temperature and daylight for base plates and individual impression trays for full or partial dentures in the dental laboratory. The material is **not suitable or recommended to be used in the patient's mouth (intraoral) for a longer time**, and only when fully polymerized. The individual impression tray is to be formed on a die model and polymerized by an appropriate fluorescent lamp during 3 – 5 minutes (wave-length: 350-400 nm). A longer polymerization time improves the physical material properties.

Any claims on warranty or liability will not be accepted if the product is not used according to determinations !

3. Caution:

Ensure that unqualified persons (e.g. children) do not come in contact with the product. Only skilled authorized personnel is entitled to work with this product. Read operating instructions before use and in case of obscurities ask the manufacturer or distributor.

4. Operating instructions:

Open the box and take out the number of light curing sheets required. Remove separating material from the sheets and close the box with the remaining material. Mould the light curing material onto the the die model wich has been isolated previously (by wax or isolation die against plastic) and cut the periphery. Then polymerize for 3-5 minutes. Take the rest of material to mould the handle which can be adapted – step by step - before or after polymerizing. The inhibition film caused by oxygene can be removed by a tissue moistened with alcohol.

5. Disposal: According to current local regulations of responsible disposal authorities.

6. Durability: 24 months after opening the original packaging.

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