www.heilblau.com 921705 IN 07/05

 **artex<sup>®</sup>cn**

Artikel-Nr. 217310

**Registrieren | Artikulieren**  
Registration | Articulation  
Enregistrement de l'occlusion  
Montage sur articulateur  
Registrazione | Messa in  
articolatore  
Registro | Articulación

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**User Manual**

**Mode d'emploi**

**Istruzioni d'uso**

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Deutsch	05 - 05
English	06 - 06
Française	07 - 07
Italiano	08 - 08
Español	09 - 09

Made in the European Union

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## MANUAL ARTEX®CN

### OPEN CENTRIC

Swivel (1) locking bar upwards with the thumb; the top part of the Artex can be guided backwards freely. The locking bar limits travel in a posterior direction and prevents the top part of the Artex falling out. Bennett angle can optionally be set from 0° to 20°.

Centric open (2), condyle ball released.

### CLOSE CENTRIC

Press locking (3) bar down; the plastic tongue centres the condyle ball in the lowest position of the condyle path. Only rotation is possible with the top part.

Centric closed (4): the condyle ball is fixed in position by the plastic tongue.

### SETTING BENNETT ANGLE

Open the Bennett (5) setting screw, turn condyle to desired Bennett angle and fix position with the setting screw.

Fig.: Bennett angle = 0°.

Bennett angle (6) can optionally be set from 0° to 20°.

Fig.: Bennett angle = 20°.

#### Accessories/single components:

Order No.	
217330	Base Plates blue, pair
217331	Base Plates blue, 50/Pkg.
215250	Retention Disk

## ! WARNINGS AND SAFETY ASPECTS:

### Anatomical anterior guidance

Angle of protrusion track = 40°;  
on lateral movements the angle varies between 30 and 40°, according to the predetermined Bennett angle

### How to check the centric adjustment of Artex CN and Artex CT

If the upper part of the articulator has lateral play, this deviation to the side can be heard and felt by a "click". In this case, it is necessary to readjust the centric.

Then, please proceed as follows:

- Set the Bennett angle to 10°.
- Close the centric lock.
- Unscrew the thread pins at the bottom of the rear part using an allen wrench.
- In the anatomically shaped anterior guidance plate, the anterior guidance pin is centrally adjusted at zero.
- Hold the upper articulator part in this position. At the same time, pull the left semi-axis outwardly until reaching the stop, and use the allen wrench to tighten the thread pin. Now proceed in the same way with the right semi-axis.

### General information

Articulators are precision instruments. Therefore, it is important to handle them with care.

Please read the following notes:

- Periodically check if the axis-joint system and the condylar paths are clean in order to maintain the smooth mobility of the condyle in the axis and thus maintain the functionality of the articulator.
- Remove dirt in this area with a steamcleaner, then blow dry.
- Do not place the articulator in a boil-out unit, a pressure pot or an ultrasonic bath.
- For cleaning and care, use only cleaning agents that are suitable for light metals. Regarding the use of the cleaning agent, please follow the instructions of the according manufacturer. Do not use cleaning agents with a pH <5 or >8 since those might damage the anodized surfaces.

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