**FINOGUM PREMIUM/
FINOGUM IMPLANT/
FINOGUM SCAN**
Zahnfleischmaskenmaterial

REF 19093/19094/19098

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REF 19093/19094/19098

ENG Instructions for use

Please read carefully before use
Keep for referral purposes

Only to be used by qualified
personnel!

1. Fields of application

Addition curing, gingiva-coloured vinyl polysiloxanes.

FINOGUM PREMIUM with low final hardness: for the fabrication of flexible gingival masks in crown and bridge techniques.

FINOGUM IMPLANT with high final hardness: specifically for implant techniques.

FINOGUM SCAN with medium end hardness: scannable and for universal use with all conventional CAD/CAM systems.

2. Safety precautions

For safety-relevant data please see the safety data sheet.

3. Processing instructions

Note:

When using a new cartridge, discard the initial volume (pea-size) to ensure uniform dispensing of the two components. Never pass the material through a hardened or hard mixing cannula, this can lead to leaks in the cartridge system and the discharge of unmixed, non-curing individual components. Check the cartridge for uniform dispensing prior to every application. Latex gloves, hand creams or cleansing agents can negatively affect the setting behaviour of the gingival mask.

3.1 Indirect method
FINOGUM PREMIUM/FINOGUM SCAN

3.1.1 Preparation

Place a matrix of hard silicone modelling material (e.g. FINOPASTE) on the solid master model where the gingival mask is to be prepared. Trim the matrix after curing such that the stump tips protrude by approx. 1 mm. Saw the model segment and preparation margin as usual. Check the exact fit of the matrix and insulate the silicone putty material with the FINOGUM separating agent included in the scope of delivery to ensure a perfect surface of the gingival mask and reliable separation of the two silicone materials.

3.1.2 Application

Apply FINOGUM PREMIUM/FINOGUM SCAN directly into the matrix and onto the master model. Place the matrix back onto the model and remove the excess escaping from the stump tips.

3.1.3 Curing

Ideally in a pressure pot. Remove the matrix after 6 min 30 s.

3.1.4 Post-processing

Post-processing of the gingival mask with scalpel, scissors or rotating instruments.

3.2 Direct method
FINOGUM IMPLANT/FINOGUM SCAN

3.2.1 Preparation

To achieve a perfect surface, apply the separating agent included in the scope of delivery to the contact surfaces of the impression material.

Note:

Use of the FINOGUM IMPLANT/FINOGUM SCAN separating agent is absolutely essential when using the gingival mask material in conjunction with polyether impressions. The polyether impression must be at least 3 h old before applying FINOGUM IMPLANT/FINOGUM SCAN.

3.2.2 Application

Apply FINOGUM IMPLANT/FINOGUM SCAN directly from the cartridge into the impression in the area of the impression posts.

3.2.3 Curing

Curing time 6 min 30 s. Then pour plaster into the impression as usual.

3.2.4 Post-processing

Post-processing of the gingival mask with scalpel, scissors or rotating instruments.

4. Storage

P411+P235 Store at temperatures not exceeding 15-28 °C. Keep cool.

5. Physical data

At room temperature of 23 °C and normal relative humidity of 50 %. Lower temperatures prolong the processing time, higher temperatures shorten the processing times.

FINOGUM PREMIUM

Processing time 1 min 30 s; tensile strength approx. 4.5 N/mm²; final hardness 65 Shore A; elongation at break 175 %; demoulding after 6 min 30 s.

FINOGUM IMPLANT

Processing time 1 min 30 s; tensile strength approx. 3.9 N/mm²; final hardness 75 Shore A; elongation at break 55 %; demoulding after 6 min 30 s

FINOGUM SCAN

Processing time 1 min 30 s; tensile strength approx. 4.5 N/mm²; final hardness 65 Shore A; elongation at break 175 %; demoulding after 6 min 30 s.

5.1 Composition

Vinyl-containing polydimethyl siloxanes, polymethyl hydrogen siloxanes, silicic acid, fillers, pigments, Pt-catalyst.

6. Delivery forms

FINOGUM PREMIUM

2 FINOGUM PREMIUM double cartridges each 50 ml; 12 mixing tips, light blue;

1 FINOGUM PREMIUM separating agent 10 ml.

1 Set 19093

FINOGUM IMPLANT

2 FINOGUM IMPLANT double cartridges each 50 ml; 12 mixing tips, light blue;

1 FINOGUM IMPLANT separating agent 10 ml.

1 Set 19094

FINOGUM SCAN

2 FINOGUM SCAN double cartridges each 50 ml; 12 mixing tips, light blue; 1 FINOGUM SCAN separating agent 10 ml.

1 Set 19098

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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