3M ESPE**ESPE™ Sil****Silane Coupling Agent**

- (D) **Silanolösung**
- (F) **Solution de silane**
- (I) **Soluzione silanizzante**
- (E) **Solución de silano**
- (P) **Solução de silano**
- (NL) **Silaanvloeistof**
- (GR) **Διάλυμα σιλανίου**
- (S) **Silanlösning**
- (FIN) **Silaaniliuos**
- (DK) **Silanoplösnings**
- (N) **Silanlösning**

Instructions for Use
Gebrauchsinfo
Mode d'emploi
Informazioni per l'uso
Informaciones de uso
Instruções de uso

3M ESPE AG
Dental Products
D-82229 Seefeld - Germany

3M ESPE
Dental Products
St. Paul, MN 55144-1000

3M, ESPE, CoJet, Rocatec, Sinfony and Visio are trademarks
of 3M or 3M ESPE AG.

© 3M 2001. All rights reserved.

Gebruiksinfo
Πληροφορίες Χρήσης
Bruksanvisning
Käyttöohjeet
Brugsanvisning
Bruksveiledning

768051/01

ENGLISH**Product Description**

Optimal adhesion between metal – composite, ceramic – composite or composite – composite is provided by a chemical bond. This bond can be established by silicatization – with Rocatec™ Plus or Soft or Cojet™ Sand – and subsequent silanization with ESPE™ Sil Silane coupling agent.

Rocatec Plus and Soft and ESPE Sil are part of the Rocatec™ System; Cojet Sand and ESPE Sil are components of the Cojet™ System, all products manufactured by 3M ESPE.

For details on all products mentioned below, please refer to the corresponding Instructions for Use.

Instructions for Use should not be discarded for the duration of product use.

Areas of Application

- Silanization of silicatized metal frameworks prior to veneering
- Silanization of silicatized ceramic, composite or metal restorations prior to cementation
- Silanization of silicatized defects in ceramic or composite veneers with or without exposed metal parts, all-ceramic or all-composite restorations or composite fillings

Preparation

In intraoral applications place a rubber dam.

Use Rocatec Plus/Soft or Cojet Sand to silanize the surface to be silanized.

Application

For each use, freshly dose ESPE Sil in a clean, grease-free dish.

- Immediately after use reseal the bottle containing the volatile solution.

Soak the brush in ESPE Sil.

- Use only the designated brush for application of the silane coupling agent.

Wet the Rocatec Plus/Soft- or Cojet Sand-silicatized surface(s) just once.

Allow the volatile silane solution to **dry for 30 sec in intraoral applications or for 5 min if applied extraorally**.

- The material must be allowed to dry for the specified times in order to allow all volatile solution components to evaporate, and prevent bonding defects!

- The silanized surfaces must be completely dry before applying opaquer or Visio™ Bond, manufactured by 3M ESPE, and composite.

- Avoid contaminating the silanized surfaces.

Immediately after drying, proceed as required in the application at hand.

Comply with the Instructions for Use of the individual products:

- Metal frame veneering: apply Sinfony™ Opaquer, manufactured by 3M ESPE.

- Adhesive Cementation: proceed in accordance with the Instructions for Use of the adhesive cement.

- Intraoral repair of ceramic or composite restorations with exposed metal parts: apply Sinfony Opaquer to the exposed metal surface and light-cure with a halogen light unit. Apply Visio Bond to the ceramic or composite and cure with light.

- Intraoral repair of ceramic or composite restorations with no exposed metal parts: apply Visio Bond to the ceramic or composite and cure with light.

Then finish the repair with composite complying with the respective manufacturer's instructions for Use.

Incompatibilities

In susceptible individuals, sensitization to the product cannot be excluded. Use of the product should be discontinued and the product completely removed, if allergic reactions are observed.

Storage and Shelf-life

Do not store the product above 25°C/77°F.

Do not use after the expiration date.

Customer Information

For questions or comments in U.S.A. or Canada please call toll-free 1-800-634-2249.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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