 **FINOSCAN**
Spray 2

REF 55552

DEU **Gebrauchsinformation**

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REF 55552

ENG **User Information**

Please read the entire document before using the Marking.

1. Intended purpose
Titanium dioxide-free, micro-pulverising scan spray.

2. Indication
For improving the optical properties of prepared models when recording by camera or scanner with the CAD/CAM technique.
For extraoral use, exclusively for indirect application in the laboratory.
Universally applicable for all CAD/CAM systems.

2.1 Contraindication
Do not use if incompatibilities or allergies to constituents have been confirmed.

3. Fields of application
For processing in the dental technology according to the intended use.
Only for extraoral use!

4. Safety precautions

  Danger

H222 Extremely flammable aerosol.
H229 Pressurized container: may burst if heated.
H411 Toxic to aquatic life with long lasting effects.

4.1 Prevention
P102 Keep out of reach of children.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P211 Do not spray on an open flame or other ignition source.
P251 Do not pierce or burn, even after use.
P260 Do not breathe aerosol.
P271 Use only outdoors or in a well-ventilated area.
P273 Avoid release to the environment.
Ensure good ventilation and venting of the work areas.

5. Processing instructions

 Shake well before use!
Place the spray nozzle on the spray can.
Clean and dry the surface to be sprayed well.
Hold the can vertically and spray onto the dry, clean surface from a distance of 5.0 - 8.0 cm.
Fully press down the spray head!
Only start recording after drying. After recording, the surface can be easily cleaned with water or air.

6. Storage
P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.

7. Physical data
Boiling range -48 °C; Ignition temperature 365 °C; Vapour pressure 2.1 hPa; Density 0.64 g/cm³; max. VOC content >95 %; Grain size 1 to 50 µm.

7.1 Composition
titanium dioxide-free.

8. Delivery forms
Packaging units for Marking are:

- Spray 2
- User Information

FINOSCAN Spray 2 with 2 spray heads net 50 ml (gross 75 ml) 55552

Spray 2 are shown in the illustration:

9. Guarantee
Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

10. Exclusion of liability
All rights reserved. This document may not be reproduced, transmitted, transferred, distributed, modified, summarised or translated into any other language graphically, electronically or mechanically, including but not limited to photocopying, recording, information storage and processing in systems, without the prior written consent of FINO GmbH. Copies of the system software are not permitted.

11. Disposal
P501 Dispose of contents/container in accordance with local/regional/national/ international regulations.

11.1 EWC waste code product
160504 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; gases in pressure containers (including halons) containing hazardous substances; hazardous waste

12. Symbols explanation

12.1 Manufacturing

REF Item numbers
LOT Lot code

Shelf life
 non-sterile
 Unit in the packaging
 Manufacturer

12.2 Application
General danger symbol warns of the risk of serious injury in conjunction with the warning words CAUTION, WARNING or DANGER

Notes
Important information
Please read carefully before use
Keep for referral purposes
Do not reuse

12.3 Storage
Temperature limit
Keep away from direct sunlight
Protect against moisture
top - store upright
Caution fragile

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