SRN: DE-MF-000021120

FINOHIT ZR HT/ FINOHIT ZR ML

Basis UDI-DI: ++EWPDCZRMV

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 206458-0

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • Germany • +49 9708 909-420

	Gebrauchsanweisung	2
	Instructions for use	5
	Manuel d'utilisation	8
	Istruzioni per l'uso	11
	Instrucciones de uso	14
	Instrukcje użytkowania	17
	Instruktioner för användning	20
	Instruktioner til brug	23
	Bruksanvisning	26
	Käyttöohjeet	29
	Gebruiksaanwijzing	32
	Οδηγίες χρήσης	35
	Upute za uporabu	38
	Instrucțiuni de utilizare	41
	Naudojimo instrukcijos	44
	Návod k použití	47
	Návod na použitie	50
	Használati utasítás	53

SRN: DE-MF-000021120

FINOHIT ZR HT/ FINOHIT ZR ML

Basis UDI-DI: ++EWPDCZRMV

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 206458-0

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • Germany • +49 9708 909-420

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Produkt/ Produkttyp

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Instructions for use

Product/ Product type

Pre-sintered zirconia blanks for the fabrication of individual dental restorations

Material type

ZrO₂ (yttria-stabilised tetragonal zirconia)/ ceramic type 2, Class 5,
Class IIa medical device

User

Trained users who fabricate individual dental restorations from zirconia.

Indication/ purpose

The devices are medical devices intended for the manufacture of dental prostheses for long-term use, which are partially introduced into the human body by means of a clinical intervention, but in this case are inserted into the teeth or fixed to the tooth substance with suitable fixing material and are thus to be classified as class IIa.

Basic frameworks for crowns, bridges, etc., as well as fully anatomical restorations for dental technology can be milled. The restorations can be veneered with ceramic and/or glazed.

Crowns and bridges up to 16 units with max. 2 pontics between 2 crowns in the posterior region, with max. 4 pontics between 2 crowns in the anterior region, veneers, inlays, onlays, primary telescopes, bar constructions

Contraindication

Do not use in case of proven hypersensitivity to one or more ingredients.
Do not use if there is insufficient space.

Veneering ceramics/ glaze

All ZrO₂ veneering ceramics/ glazes

Material specific properties

Chemical composition	FINOHIT ZR HT A1 - D4	FINOHIT ZR ML A1 - D4
ZrO ₂ + Y ₂ O ₃	≥98,7	≥98,7
Y ₂ O ₃	6,66 (± 0,25)	6,66 (± 0,25)
Al ₂ O ₃	0,05 (± 0,01)	0,05 (± 0,01)
Fe ₂ O ₃	0,08 (± 0,07)	0,08 (± 0,07)
ER ₃ O ₃	≤ 0,56 %	≤ 0,56 %
CO ₃ O ₄	≤ 0,01 %	≤ 0,01 %
Other oxides	≤ 0,02 %	≤ 0,02 %

Physical properties	FINOHIT ZR HT A1 - D4	FINOHIT ZR ML A1 - D4
Density g/cm ³	6,07 - 6,33	6,07 - 6,33
Bending strength	~1100 MPa	~1100 MPa
Coefficient of thermal expansion (25°C -500°C)	10.5 +/- 0.5 x 10 /K ⁻⁶	10.5 +/- 0.5 x 10 /K ⁻⁶

Processing

The desired restoration is milled from the selected blank:

SRN: DE-MF-000021120

FINOHIT ZR HT/ FINOHIT ZR ML

Basis UDI-DI: ++EWPDCZRMV

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 206458-0

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • Germany • +49 9708 909-420

The wall thickness for anterior teeth must not be less than 0.5mm and for posterior teeth 0.6mm. The junction between two splinted anterior crowns must not be less than 8mm². The connection between an anterior crown and a pontic, two splinted posterior crowns or two anterior pontics must not be less than 9mm². The connection between a posterior crown and a pontic or two posterior pontics must not be less than 12mm². The cross-section of a bar construction must not be less than 15mm .²

The dimensions refer to the dimensions after sintering for all product variants.

The enlargement or shrinkage factor is indicated on the blank depending on the milling system used.

The unsintered restorations are still finished if necessary. In order to obtain a perfect surface and maximum translucency, we recommend cleaning the restoration in the ultrasonic cleaning unit "White-Sonic" with distilled water. Please use only our recommended ultrasonic cleaning device. Devices from other companies could be too "strong" and damage the restoration. Swivel the restoration with plastic tweezers in the distilled water of the ultrasonic unit for 5 to 10 seconds until no more "dust cloud" comes off the restoration. Remove excess water from the restoration by blowing it off with oil-free air or by blotting it dry with cellulose or cotton swabs. Dry the restoration under an infrared lamp or using an oven. Before sintering, the restoration must be completely dry. Make sure that the water never boils during drying, as this may cause cracks in the zirconia.

The finished, completely dry restoration is now sintered.

Sintering

	Normal Program	Slow Program	Translucency Program	Speed Program
Heating rate	10°C per minute to 950°C	5°C per minute to 950°C	5°C per minute to 950°C	50°C per minute to 1100°C
Holding time	none	none	none	none
Heating rate	6°C per minute to 1500°C	2°C per minute to 1500°C	2°C per minute to final temperature	20°C per minute to 1500°C
Holding time	at final temperature 90 minutes	at final temperature 120 minutes	at final temperature 120 minutes	at final temperature 30 minutes
Final temperature	1500°C	1500°C	1500°C - 1630°C	1500°C
Cooling down	unregulated in the closed oven	unregulated in the closed oven	unregulated in the closed oven	unregulated in the closed oven

Post-processing

After sintering, the restoration is adjusted to the working model, if necessary, by wet grinding with diamond-coated grinding instruments. Sintered diamonds, corundum stones or carbide burs should not be used. Overheating must be avoided.

Firing the ceramic/glaze

All commercially available firing-on ceramics/glazes for zirconium dioxide can be used. Please observe the working instructions of the respective manufacturer.

Firing and cooling according to the firing chart of the manufacturer of the ceramic mass/glaze used.

Safety instructions

Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and clothing and switch on extraction equipment. The dust generated by this product during processing may cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/ contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after use. Remove contaminated clothing and protective equipment before entering areas where food will be eaten. Keep away from sources of ignition. Do not smoke .

SRN: DE-MF-000021120

FINOHIT ZR HT/ FINOHIT ZR ML

Basis UDI-DI: ++EWPDCZRMV

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 206458-0

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • Germany • +49 9708 909-420

Storage

Store in a dry place. Protect from moisture. Store in the original packaging.

Disposal

Dispose of product and packaging according to local/ regional/ national/ international regulations. Do not dispose of together with household waste. Do not allow to enter waters, groundwater or sewage system.

Note

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

The summary report on safety and clinical performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European Database for Medicinal Products (Eudamed) under the associated base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

Explanation of the markings on the packaging

Confirmation: The product complies with the applicable European directives

Catalogue number/ Item number

Batch designation

Clear product identification

Number of products in the package

Please follow the instructions for use

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